

Introduction

NeuroOptics offers a wide selection of infrared imaging pupillometers to meet your research needs. The questions to answer to narrow down your options are as follows:

- Human or animal?
- Monocular (hand-held) or binocular?
- Pupil diameter only or also the pupillary light reaction?
- Need to download data?

VIP®-300 Pupillometer

The NeuroOptics VIP-300 Pupillometer is a portable, battery-operated, hand-held device which very accurately measures pupil size under various light conditions.




Ordering Information

NeuroOptics® VIP®-300	Part Number
System Includes: VIP®-300 Pupillometer, Charging Station & Power Supply, Eye Cups (2), Carrying Case	VIP3-KIT-01
Optional Accessories	Part Number
Printer Kit	NEUR-PRTS445-BT
Barcode Scanner	BCS-CC-04-1D

PLR[®]-3000 Pupillometer

The PLR-3000 Pupillometer is a portable, battery operated, hand-held device which quickly and accurately measures pupil size and the pupillary light reflex. Used by researchers at leading universities, pharmaceutical companies, and contract research organizations worldwide, the PLR-3000 is an important research tool.




Ordering Information

NeurOptics [®] PLR [®] -3000	Part Number
System Includes: PLR [®] -3000 Pupillometer, Charging Station & Power Supply, Eye Cups (2), Carrying Case	PLR3-KIT-01
Optional Accessories	Part Number
Printer Kit	NEUR-PRTS445-BT
Barcode Scanner	BCS-CC-04-1D

DP-2000 Pupillometer

The DP-2000 Pupillometer is a binocular, multi-chromatic (white, green, blue and red) system which stimulates direct, consensual, or both eyes simultaneously. The DP-2000 is a very sophisticated research tool for many applications.




Ordering Information

NeurOptics [®] DP-2000 Pupillometer and Laptop	Part Number
System Includes: DP-2000 Pupillometer and Laptop	DP-2000-KIT-01

A-2000 Small Animal Pupillometer

The A-2000 Pupillometer is a binocular, multi-chromatic (white, green, blue and red) system which stimulates direct, consensual, or both eyes simultaneously. The device is designed to be used with a variety of animals including rodents, rabbits, and other animals with similar pupils.



Ordering Information



NeuroOptics® A-2000 Small Animal Pupillometer

Part Number

System Includes: A-2000 Pupillometer and Laptop

A-2000-KIT-01

Which NeurOptics Pupillometer is Right for Me?



NEUR OPTICS

Pupillary Evaluation...Down to a Science™

9223 Research Drive
Irvine, CA 92618 | USA
p: 949.250.9792
Toll Free North America: 866.99.PUPIL
info@NeurOptics.com
NeurOptics.com